

Thank you!

For choosing our Gymna Mobile Fit.
We have developed this trolley with great care
and attention. We hope it will make your work
more enjoyable every day.

Thanks again for your support!

gymna



Get to know your Gymna Mobile Fit

Detachable vaco cup/cable holder

For tidy storage.

Inlay top tray

For safe placement of Gymna Ultrasound Compact or Acure devices.
Removable for easy cleaning.

Detachable drawer

Allowing easy cleaning when needed.
Inner dimensions: 225x160x75mm.

Sturdy central column

Allowing power cable concealment and protection.

A foot attachment

Allows repositioning of the Mobile Fit by using the foot.
An advantage if the physical therapist has gel or lotion on their hands.

Two handles

For maneuvering within the practice.

Detachable bottle holder

The bottle may also be stored upside down.
The holder can be mounted on either left or right side.

Detachable cable organizer

Holder for tidy storage.
Handy when moving the Mobile Fit.

Cable opening

For power cable management/guidance.

Four swivel wheels

For maneuvering freely and easily.
Lockable front wheels to hold the Gymna Mobile Fit in position.



Discover more
at gymna.com